

# EyeLink Models

---

The different rooms have slightly different models of the EyeLink devices. Here is an overview. Make sure to consider this carefully!

## Space B449 & EyeLink Cubicles

The rooms in B449 (a, b, c, and d) and the EyeLink cubicles (A437m/n & A437p/o) make use of:

### EyeLink 1000 Plus Desktop Mount

Stabilized: Up to 2000 Hz monocular and binocular. Remote: Up to 1000 Hz monocular and binocular.

## Space B447

The rooms in B447 (a, b, and c) and the EyeLink cubicles (A437m/n & A437p/o) make use of:

### EyeLink 1000 Plus Desktop Mount


Stabilized: Up to 2000 Hz monocular and up to 1000 Hz binocular. Remote: Up to 500 Hz monocular and binocular.

## Tower EyeLink System

There is also the option to use the **Eyelink 1000 Tower Mount** system (for example for experiments including touch screen). It can go up to 1000 Hz with one eye.

Please send an email to the Lab Coordinator in advance, as this system is not a default set-up. It will need to be set-up before use.

The EyeLink 1000 Tower Mount camera:



---

Revision #4

Created 2025-11-03 15:02:56 UTC by Admin

Updated 2026-06-08 13:43:25 UTC by Admin